

### Midmark 530, 535 Spec Sheet

<u>Description</u>	<u>Dimension</u>
Dimensions (Model 530) – length	84” (213.4 cm)
– width	29” (73.7 cm)
– height (Adjustable)	25.75” to 35.75 in (65.4 cm to 90.8 cm)
Mattress	78” x 25.5” (198.1 cm x 64.5 cm)
Dimensions (Model 535) – length	84” (213.4 cm)
– width	33” (83.8 cm)
– height (Adjustable)	25.75” to 35.75” (65.4 cm to 90.8 cm)
Mattress	78” x 29.5” (198.1 cm x 74.9 cm)
Siderails	12” (30.5 cm) high above litter top x 63.75” (161.9 cm) length
Casters (Central locking of the brakes and steering may be accomplished from all four corners of the stretcher.)	8” (20.3 cm)
Fowler	0° to 90°, pneumatic assisted
Trendelenburg	12°, pneumatic assisted
Reverse Trendelenburg	-7°, pneumatic assisted
Oil used in pump system (530/535 Models)	Light grade mineral oil [clear], Dexron II transmission fluid [red], or power steering fluid [light amber] (Use the same type of oil presently used in the stretcher pump – Do not mix oils.)
Maximum patient load	350 lbs (158.7 kgs)
Shipping Carton (All Models)	88” x 33” x 38” (223.5 cm x 83.8 cm x 96.5 cm)
Weight (Model 530) – normal	266 lb (120.6 kg)
Weight (Model 530) – with shipping carton	360 lb (163.3 kg)
Weight (Model 535) – normal	268 lb (121.6 kg)
Weight (Model 535) – with shipping carton	370 lb (167.8 kg)